IMPROSCAR PATCH™

IMPROSCAR PATCH™

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IMPROSCAR PATCH™

Silicone sheets for treatment of scars

Water resistant

7 days continuous wear time

EN

Directions for use: See illustrations below.

Clean your hands, the scar and surrounding skin using soap and water. Following this, clean the areas with the provided Improscar Wipe™. Once dry, apply Improscar Patch™ on the scar ensuring not to touch the adhesive side which could otherwise lead to a shorter wear time. The patch should cover the whole scar and extend approximately one centimeter beyond the scar's edges. Improscar Patch™ is designed to remain securely adhered to the skin for an average duration of one week and can be worn in the shower. Use Improscar Patch™ for a minimum of 3 months and continue for as long as the scar responds to the treatment.

When changing the silicone sheet, allow the skin to rest for a few hours until any redness disappears before applying a new Improscar Patch™.

For longer scars: Use two or more Improscar Patch™ to cover the whole scar and let the patches overlap by approximately 1 cm.

Warnings: Avoid contact with the eyes and mucous membranes.

Do not use on damaged or broken skin (open wounds).

If rash occurs, stop use and consult a doctor.

Storage conditions: Store between 10°C (32°F) - 30°C (86°F), keep away from direct sunlight.

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